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
CLONOS is a NEOORTHO's trademark that uses french company technology from Biomatlante to manufacture the Resorbable Bone Substitute.

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
KEY FEATURES	KEY BENEFITS
Osteoconductive	Provides a scaffold for new bone growth.
Molecular mixture of: 60% HA and 40% TCP	HA alone resorbs too slowly while TCP resorbs too fast. A bi-phasic of HA and TCP allows the resorption rate to be similar to that human bone.
70% porosity, interconnected network of macropores and micropores	Porosity, similar to that of cancellous bone, allows the bone cells colonization and biological fluid uniformly inside the matrix.
Microporosity (<10 microns)	For ionic exchange: TCP dissolution and bone crystal precipitation. Newly bioactive interface with bone cells.
Macroporosity (>10 microns)	Allows deep invasion of bone cells into the matrix.
Technology with more than 30 years of clinical experience	Host bone formation is well demonstrated time after time.
Safe	5 years shelf life - Fully synthetic.
Convenient	Available in different granules for all indications.

AVAILABLE:	GRANULOMETRY (mm)	QUANTITY	CODE
VIAL	0,5-1,0	0,5cc/0,35g	9.001.901
		1,0cc/0,7g	9.001.903
	1,0 - 2,0	1,0cc/0,6g	9.001.116
		2,0cc/1,2g	9.001.118
SERINGA	0,5-1,0	5,0cc/2,5g	9.001.222
		10,0cc/5,0g	9.001.232
	2,0 - 3,0	0,5cc/0,35g	9.001.900
		1,0cc/0,7g	9.001.904
2,0 - 3,0	5,0cc/2,5g	9.001.223	
	10,0cc/5,0g	9.001.231	

DIRECTIONS

- 1 - Connect the syringe holding the bone marrow to Clonos syringe
  - 2 - Inject the bone marrow previously aspirated
- Clonos is ready to use

ABOUT CLONOS

Clonos is a bioactive calcium phosphate ceramic composed of hydroxyapatite (HA) and Tricalcium Phosphate (TCP). Clonos is indicated to substitute for bone graft in nonloadbearing clinical applications and is supplied in both block and granular form.

- PERMEABLE**  
3D interconnected scaffold helping cells spread  
1/3 Micropores < 10µm  
2/3 Macropores 300-600µm
- BIOACTIVE**  
Biphasic Calcium Phosphate  
Short term dissolution and precipitation process into micropores.

- SAFE**  
100% Synthetic  
No disease transmission  
30 years of clinical background.
- RESORBABLE**  
Real chemical synthesis: Homogenous distribution of 60% HA (Hydroxyapatite) and 40% TCP (Tri Calcium Phosphate) Thin crystalline structure.
- OSTEOCONDUTIVO**  
Macropores guides bone cells into the depth of the Clonos implant.
- BONE REGENERATION**  
Simulates the osteoblastic response: proliferation and differentiation 100% regeneration into cortical or cancellous bone.